Titan[®] Spring Deck Bunk Bed

SPECIFICATIONS

DETAILS ·

*Minimums apply

Options

Assembly

Color



Steel Color (Spring Deck): Matte Ebony

Steel Color (End Frames): Khaki, Ebony

Optional Steel Colors: Cream, London Fog,

Pewter, Shark, Slate, Irish Moss, Dusk, Mink

*Due to differences in monitors, printers and materials,

actual product colors may vary. Color chips are available

Standard), Ladder, Safety rail, Step Ladder,

Attach wedge lock brackets to end frames

with furnished hex socket button head

members of end frame). Select desired

height for bed deck prior to assembly.

Secure spring deck bolts to wedge lock slots. Optional accessories attach to bed (with hardware supplied through purchase of

fasteners (threaded inserts in vertical

Headboard & footboard (Wild Cherry

for more accurate depiction of color.

Bolt-down Kit, Glide Kit

Model

TNT1351 Weight 172 lbs

Dimensions Length 80.6" 38.5" Width 72" Height **Recommended Mattress**

36"x75"

*Intended for indoor use only

USA

Fabrication & Material

Cut, formed, welded steel, powder coated steel.

Bed deck: 2" x ¹/₅" x ⁵/₃₂" rail angle steel, sinuous 9 gauge wire connected by 13 gauge links with 12 gauge 16 turn helical springs.

End frame legs: 16 gauge

End frame cross-members: 16 gauge

End Caps: Injection molded high-impact plastic

Optional Headboard: 45 lb density M3 grade particle board with durable highpressure laminate surface and 2mm PVC edge (Wild Cherry Standard). Made in the U.S.A.

Installation

None required for standard model

Warranty

10 year limited replacement warranty

Maintenance

Norix recommends the use of soap and water for both the regular cleaning maintenance and deep-scrub cleaning of powder coated steel.

Optional Decorative Laminate: Norix recommends the use of Simple Green® or an equivalent cleaning product.

TESTING

accessories).

Environmental Test Standards

• GREENGUARD and GREENGUARD Gold Certified by UL Environment. Product certified for low chemical emissions to UL 2818. UL.COM/GG.

• Compliant with Healthier Hospitals guidelines

Performance Test Standards

Tested to 800 lb static load







NORIX.COM